

Perma Hydra

Our non-invasive treatment improves skin health, addresses individual concerns and creates a glow like no other. Unlike traditional facials, it's performed via device rather than by hand, for deeper penetration and longer lasting results.



TECHNICAL SPECIFICATIONS:

Voltage	100-120V/ 60Hz 200-240V/50Hz
Display	10.4 inches
Power	500W
RF Frequency	0.2-0.5MHz
Heat&Cool	Heat: 36-46 Cool: 0-6
Ultrasound	1Mhz
Aqua Peel Probes	Soft: 3 pcs Hard: 8 pcs
Aqua Peel	Vacuum intensity: 4-5kg Vacuum Pressure: 96Kpa Water flow: ≤50L
Solution Bottle	A: Cleaning B: Whitening and moisturizing C: Deep layer cleaning D: Pipe cleaning
CLEANSING SHOVEL Frequency	26-30MHz
EMS NANO Mesotherapy Size	4X4mmX250Um
MICRO LIFT	Frequency: 200Hz Current: 50-800Ua
Skin Analyzer	28 MILLION HD MEGA PIXELS

Perma Hydra

Get the Best Skin
of Your Life



Perma Hydra

Meet all your skin care needs!

- Fine lines and wrinkles
- Elasticity and firmness
- Even tone and vibrancy
- Skin texture
- Brown spots
- Oily and congested skin
- Enlarged pores

4 STEPS. 30 MINUTES.

Get the Best Skin of Your Life

It delivers instant and long-lasting results that you can see and feel.

STEP 1 ANALYSER

200x magnifying lens for deep and clear checking of skin conditions.

STEP 2 Cleanse + Peel

Uncover a new layer of skin with gentle exfoliation and relaxing resurfacing.

STEP 3 Extract + Hydrate

Remove debris from pores with painless suction.
Nourish with intense moisturizers that quench skin.

STEP 4 Fuse + Protect

Saturate the skin's surface with antioxidants and peptides to maximize your glow.

MULTIFUNCTIONAL VERSATILE APPLICATORS

1.AQUA PEEL



Use of a vacuum to open pores for a deep cleanse and easy extraction of impurities and dead skin cells. It also allows the active ingredients in our boosters to penetrate deeper, for the most effective treatment possible.

2.OXYGEN JET



High-pressure oxygen injection with small particles and strong penetration. Penetrates deep into the basal layer, enhancing water retention, absorbing more easily by narrowing pores.

3.EMS FACE



Nanocrystals promote the absorption of skin care products by opening channels for skin absorption. Therefore, when nanocrystals are used, skin care products can be more easily absorbed by the skin.

Perma
Hydra

- 1 AQUA PEEL
- . OXYGEN JET
- 2 EMS FACE
4. COLD-HEAT
3. MICRO LIFT
6. OXYGEN
7. ULTRASOUND
- 8 RF
- . CLEANSING BLADE
10. ANALYZER
- .

4.COLD-HEAT



Quickly lower the temperature below 0 degrees for skin quality, which achieves the effect of removing redness, anti-sensitivity, firming and shaping

5.MICRO LIFT



Anion and cation EMS microcurrent alternates between positive and negative, so that cations are imported and anions are exported improving facial contours.

6.OXYGEN



Create an oxygen supply environment, activate cells, and enhance absorption. Injecting the product into the deep layer of the skin and cooperating with fruit acid skin peeling products.

7.ULTRASOUND



Ultrasound—1 to 3 million vibrations per minute, opens the pores, allowing the essence to quickly penetrate into the base layer of the skin, allowing the absorption rate to reach more than 90%.

8.RF



The probe can emit radio waves that vibrate 6 million times per second. The high-frequency vibrations are the molecules in the skin that generate heat due to friction, so as to achieve the effect of firming the skin and improving wrinkles.

9.CLEANSING BLADE



Use high-frequency vibration waves to soften the stratum corneum and deeply clean the skin dirt.

10.ANALYZER



Skin Analyzer:
The best tool to make skin analysis and offer best solutions!

HOW DOES PermaHYDRA WORK ?

Nanocrystals use nano-scale contacts, which are as thin as more than 250 nanometres and as short as 0.1 mm, which can just open the stratum corneum of the skin without touching the dermis. It not only solves the problem of skin absorption, but also achieves painless and non-invasive treatment!



Natural application and absorption

In the case of natural application and absorption of skin care products, the molecules of skin care products only penetrate into the epidermis. It is difficult to penetrate into the deeper skin layer and therefore the absorption rate is very low.

Nanocrystalline nanochip (The product can penetrate evenly into the dermis.)

Nanochips only penetrates the barrier layer of the epidermis without touching blood vessels and nerves. It is painless and non-invasive, and the absorption capacity is greatly improved.

HydraPERMA function

The microelectronic heat conduction technology opens the pores, which can increase the product absorption rate by more than 40%.

Principle of nanochip introduction


- ⊙ Nanocrystals are made of monocrystalline silicon. In addition to being tough, monocrystalline silicon can also be thin enough and short enough to only break through the barrier layer of epidermal absorption, and it's far from touching blood vessels and nerves, so you can't feel pain, nor will it bleed.
- ⊙ Nanocrystals promote the absorption of skin care products by opening the skin absorption channel, so skin care products can be more absorbed by the skin when nano crystals are introduced.




Sagging skin



Dry and water shortage



Fine water wrinkles



Gloomy



DEEP CLEANSING SKIN SYSTEM

Cleaning blackheads, acne, and mites while moisturising nutrients supplying; fining pores, increasing the skin transparency and improving skin texture.

What is special about it?

Three times the power of the operating head.




The three major areas with clear division of work can achieve 3x the effect simultaneously.

ANTI-ALLERGY


SOOTHING

REPAIR

SKIN TIGHTENING



Perma
Hydra



INTELLIGENT OPERATING SYSTEM

- Versatile Functions
- Versatile Languages
- User Friendly UI
- 10.4" display
- 10 in 1 Platform



With decades of experience within the beauty sector Optilight understands the issues you face and the ultimate sense you get with a satisfied customer.

Optilight is in the unique position of being able to give you complete impartial and professional advice offering sales, service & repair solutions for the world's biggest brands of medical laser and IPL machines.

We supply consumables and even offer preventative maintenance packages so you can focus entirely on serving your customers in your salon.

- ◆ Fully trained engineers
- ◆ Nationwide service
- ◆ Gold Silver & Bronze maintenance package options
- ◆ Genuine parts
- ◆ More cost effective than using the manufacturers direct
- ◆ Honest advice



Perma™ Range



LIFE-CHANGING RESULTS

The amazing results are both instant and long-lasting.

FINE LINES



before

after


NASOLABIAL FOLDS



before

after

BROWN SPOTS



before

after

OILY & CONGESTED SKIN



before

after

ENLARGED PORES



before

after

EVEN TONE



before

after

SAGGING SKIN



before

after

DRY SKIN



before

after



0204 603 5882



optilight.co.uk



sales@optilight.co.uk



Grange Farm
Business Park,
Old Hurst Road,
Woodhurst Cambs,
PE28 3BQ